Go to traceablelive.com to create a new account.






Register Account

THEN Download TraceableLIVE mobile app via Apple App Store or Google Play.

Add Device to Account



Open TraceableLIVE app on mobile or desktop and Log in to account:



Tap the + icon to add a device.

- Permission must be given to TraceableLIVE app to access your mobile device's camera.
- Using your mobile device's camera hover over QR code on back of unit (see right).



App will automatically fill in device information.

Note: QR code only works with a mobile device.

WiFi Setup





Press/HOLD WiFi button "AP" will appear on display.

Press WiFi button again, "RP URIF" will appear on display.

When "RP rERdy" appears, proceed to next step.

Use Internet enabled device to connect to WIFI





Select Network: CC6500-0013 see example network

Example

Network: CC6500-0013



Open Browser, type 192.168.1.1 in address bar



Select local network and enter password.

Security type WPA2 recommended or contact your network admin.



CONGRATULATIONS! Your device will restart.

Still having trouble?
Please give us a call
281-482-1714 or email
support@traceable.com

TraceableLIVE® QUICK START GUIDE