


1. Eyepiece
2. Battery cover
3. Focusing Wheel
4. Light Condenser
5. Light Switch
6. Bulb
7. Transparent Cover

OPERATION

1. Place the transparent cover (7) on the object.
2. Position the hole over the area to be viewed.
3. Slide the light switch (5) to the ON position and slide down the light condenser (4) to concentrate the light on the object being viewed.
4. Look through the eye piece (1) and turn the focus wheel (3) until the object is in sharp focus.




5. When viewing is completed, turn the light switch (5) to the OFF position and retract the light condenser (4).

BATTERY REPLACEMENT

1. Remove the battery cover (2) by sliding it in the direction indicated by the arrow.

2. Remove the exhausted batteries and replace with two new AA batteries.


3. Replace the battery cover.



BULB REPLACEMENT

1. Remove the transparent cover (7) by sliding it in the direction shown in the illustration. (Note: The light condenser (4) must be retracted)

2. Unscrew the exhausted bulb and replace with a new one.
3. Replace the transparent cover.



WARRANTY OR SERVICE

For warranty or service contact:

TRACEABLE® PRODUCTS

12554 Old Galveston Rd. Suite B230
Webster, Texas 77598 USA
Ph. 281 482-1714 • Fax 281 482-9448
E-mail support@traceable.com • www.traceable.com

Traceable® Products are ISO 9001:2018 Quality-
Certified by DNV and ISO/IEC 17025:2017
accredited as a Calibration Laboratory by A2LA.

Cat. No. 3355 & 3356

©2018 Traceable® Products. 92-3356-00 Rev 1 061118

**ILLUMINATED
MICROSCOPE
INSTRUCTIONS**