| SPECIFICATIONS      |   |
|---------------------|---|
| Display:            | ½"-High LCD showing minutes and seconds   |
| Accuracy:<br>Timing | 0.01%   |
| Capacity:           | Count down – 1 second to<br>99 minutes, 59 seconds<br>Count-up – 1 second to<br>59 minutes, 59 seconds. |

Audible alarm sounds for one minute Countdown timing and

Functions: count-up (stopwatch) timing

Attachments: Clip, magnetic back, and

stand

Size:

CDECIFICATIONS

Alarm:

Weight:

21/3" x 21/4" x 1/3" 1.5 oz

the STOP/RESET button to return the display to 00.00.

15 seconds, press keys 1 and 5.)

. Press the STOP/RESET button to set

2. Press the numbers keys 0 through 9

to set desired minutes and seconds.

**Note:** To time in seconds only, press just

two of the number keys. (Example: To set

3. If the wrong time has been entered, press

(Example: To set 12 minutes 38 sec-

timer to zero. The display should read

COUNTDOWN TIMING

onds, press 1, 2, 3, 8.)

down timina.

00:00.

4. Once the desired time is on the display,

press the START button to begin count-

60 second alarm sounding time and after the alarm is turned off.

00.00.

6. During the countdown process, the timer may be paused by pressing the STOP/ RESET button one time and restarted by pressing the START button.

button or number key.

5. When the timer reaches 00:00, an alarm

will sound. The alarm will sound for 60

seconds and then turn off automatically.

To turn the alarm off manually, press any

**Note:** The display will read 00.00 during the

Note: Pressing the STOP/RESET button twice, will clear the display and return it to

count-up timing.

00.00.

Note: Timer will count up to 59 minutes, 59

seconds and then automatically reset itself to zero - the display will read 00.00 for one

second and then start counting up 3. During the count-up process, the timer

STOPWATCH (COUNT-UP) TIMING

. Press the STOP/RESET button to set the

timer to zero. The display should read

2. Press the START button once to begin

may be paused by pressing the STOP/

RESET button one time and restarted by pressing the START button. **Note:** Pressing the STOP/RESET button

twice, will clear the display and return it to 00.00.

4. To end the count-up process at any time, press the STOP/RESET button twice.

## **TIMEOUT**

The timing process may be stopped at any time during count-up or countdown. Press the STOP/RESET button one time. Timing may be resumed by pressing the START button.

**Note:** Pressing the STOP/RESET button twice, will clear the display and return it to 00.00.

## **ALL OPERATIONAL DIFFICULTIES**

If this timer does not function properly for any reason, please replace the battery with a new high quality battery (see "Battery replacement" section). Low battery power can occasionally cause any number of "apparent" operational difficulties. Replacing the battery with a new fresh battery will solve most difficulties.

## **BATTERY REPLACEMENT**

An incorrect display, no display or operational difficulties indicate that the battery should be replaced. Use a coin to open the battery cover on the back of the timer (turn the cover approximately 1/8 of a turn counter-clockwise). Equivalent battery replacements are: RAYOVAC RW42, DURACELL D357, GPA76, and EVEREADY 357. Insert the battery with the positive side facing you. Replace the battery cover. Replacement battery Cat. No. 1039.

WARRANTY, SERVICE, RECALIBRATION For warranty, service or recalibration contact:

## TRACEABLE® PRODUCTS

12554 Old Galveston Rd. Suite B230 Webster, Texas 77598 USA Ph. 281 482-1714 • Fax 281 482-9448 E-mail support@traceable.com www.traceable.com

Traceable® Products are ISO 9001:2018 Quality-Certified by DNV and ISO/IEC 17025:2017 accredited as a Calibration Laboratory by A2LA.

Traceable® is a registered trademark of Cole-Parmer.
© 2020 Traceable® Products. 92-5020-30 Rev. 4 042420

Mini-Alarm-Timer/ Stopwatch Instructions