

Go to traceablelive.com to create a new account.

1

TraceableLIVE

Email

Password

New Account

Forgot Password?


Register Account

THEN Download TraceableLIVE mobile app via Apple App Store or Google Play

Add Device to Account

2

Open the TraceableLIVE app and Log in to account:



Tap the **+** icon to add a device.

- Permission must be given to TraceableLIVE app to access your mobile device's camera.
- Using your mobile device's camera hover over QR code on back of unit (see right).
- App will automatically fill in device information.



192103821

WiFi Setup

3




Press/HOLD WiFi button "AP" will appear on display.

Press WiFi button again, "AP URIF" will appear on display.

When "AP rERdy" appears, proceed to next step.

Use Internet enabled device to connect to WIFI


4



Select Network: **CC6500-0013**
see example network

Example

Network: **CC6500-0013**



Open Browser, type
192.168.1.1 in address bar

5

CONGRATULATIONS!
Your device will restart.

Select local network and enter password.

Security type WPA2 recommended or
contact your network admin.



Still having trouble?
Please give us a call
281-482-1714 or email
support@traceable.com

TraceableLIVE®
QUICK START
GUIDE